

soga®



Operating Specification
Dual-heads Massager

The executive standard:GB4706.1-2005 GB4706.10-2008
Quality certificate:ISO9001:2008

Please hold the manual well.
Please read the manual carefully before using,especially the part of caution concerned,and keep the operation book well, please.Thank you.

Massaging Principle

Massaging is one of treatments of medicine outside .It can act on the body at specific points by all kinds of technique ,so that it can adjust system physiological and pathology condition. And then you will have a good therapy and strong body.The massaging principle mainly embodies at following points:the changes of system fuction,the ajustment of biological information and unusual correcting the postion of dissection,etc.It is main treatment means or auxiliary measures for some disease,and it is used widely.

According to the theory ,especially theory of main and collateral channels,it emphasizes that body surface of human is connected with the internal organs by main and collateral channels or points.If the internal organs are something wrong,the disease will reflect on the body surface by main and collateral.And then,we can massage and stimulate the main and collateral or points of body surface,it will make gaseous pulse move well. So,when massaging,the treatment "information" will be sent to the sick internal organs,it will make your internal organs' function stronger.

From the anatomical knowledge of modern medicine,massaging will promote the natural power of resisting disease, promote absorbing of inflammation exudation,and can dull muscle spasm and ache.

Through practice validation after thousands of years,massaging has a good treatment for inflammation,conglutination, spasm,hyperplasia,degeneration, fibrosis,blood vessels nerve of hurting parenchyma.

I. Uses and benefits

1. Ease tense muscles of back and waist.
2. Relax stiff muscles and tault nerve,especially after exercise.
3. Relieve hip,shoulder and neck's pain.

II. Function of Product

1. PTo eliminate fatigue and relax muscle because of violent action.
2. To eliminate ache because of long bad blood circulation.
3. To eliminate spasm shoulers because of stiff neck.
4. To eliminate ache of tiredness or rheumatism.
5. The spherical massage head can massage all the parts of body.
6. It can make the muscles to produce fast rhythm exercise.

III. Caution

1. Do pull out the plug after using or before cleaning it.
2. Don't use it at moist place.
3. When children or the disabled use it,pay more attention to it.
4. If the product is something wrong,do not repair it yourself or don't go on using it,you must contact with the seller ASAP.O

IV. Uses



soga®



Operating Specification
Dual-heads Massager

The executive standard:GB4706.1-2005 GB4706.10-2008
Quality certificate:ISO9001:2008

Please hold the manual well.
Please read the manual carefully before using,especially the part of caution concerned,and keep the operation book well, please.Thank you.

Massaging Principle

Massaging is one of treatments of medicine outside .It can act on the body at specific points by all kinds of technique ,so that it can adjust system physiological and pathology condition. And then you will have a good therapy and strong body.The massaging principle mainly embodies at following points:the changes of system fuction,the ajustment of biological information and unusual correcting the postion of dissection,etc.It is main treatment means or auxiliary measures for some disease,and it is used widely.

According to the theory ,especially theory of main and collateral channels,it emphasizes that body surface of human is connected with the internal organs by main and collateral channels or points.If the internal organs are something wrong,the disease will reflect on the body surface by main and collateral.And then,we can massage and stimulate the main and collateral or points of body surface,it will make gaseous pulse move well. So,when massaging,the treatment "information" will be sent to the sick internal organs,it will make your internal organs' function stronger.

From the anatomical knowledge of modern medicine,massaging will promote the natural power of resisting disease, promote absorbing of inflammation exudation,and can dull muscle spasm and ache.

Through practice validation after thousands of years,massaging has a good treatment for inflammation,conglutination, spasm,hyperplasia,degeneration, fibrosis,blood vessels nerve of hurting parenchyma.

I. Uses and benefits

1. Ease tense muscles of back and waist.
2. Relax stiff muscles and tault nerve,especially after exercise.
3. Relieve hip,shoulder and neck's pain.

II. Function of Product

1. PTo eliminate fatigue and relax muscle because of violent action.
2. To eliminate ache because of long bad blood circulation.
3. To eliminate spasm shoulers because of stiff neck.
4. To eliminate ache of tiredness or rheumatism.
5. The spherical massage head can massage all the parts of body.
6. It can make the muscles to produce fast rhythm exercise.

III. Caution

1. Do pull out the plug after using or before cleaning it.
2. Don't use it at moist place.
3. When children or the disabled use it,pay more attention to it.
4. If the product is something wrong,do not repair it yourself or don't go on using it,you must contact with the seller ASAP.O

IV. Uses




Fault	Reason	Fault Remedy
Motor not working	Plug not fixed up	Fix up the plug well.
	Switch is out of position	Check the Switch

IX. Jugegment of Fault

Power: 25W
Using time:15m/time

VII. Technological Parameter

There are around 650 points in our body from neck to waist.And every point is connected with some organs or body condition.Massaging can promote blood circulation,promote metabolism and build up the ability agaisit disease.Massager will make you have massaging delight and healthy pleasure by easy ways.



VII. Chart of body point

1. Put the plug into the socket.
2. Press "ON/OFF" button to start or close.
3. Press "VITESSE" button to increase or decrease the massage strength, four speeds to choose.
4. Press "PROGRAMME" button to choose automatic mode,three modes (RELAXANT,ENERGISANT, STIMULANT) to choose.
5. Press "HEAT" button to start or close heating function, the indicator light will be on or off. (Need to fit on heating massage head before open heating function,otherwise no heating effect.)
6. All applications will be automatically closed after 15 minutes.
7. Please pull out the plug after using.

VIIHow to use

Note: Only with heating massage head can be used for heating function



The massager hammer is equiped with three sets massage head accessories.

V. Massage Head Accessories


Fault	Reason	Fault Remedy
Motor not working	Plug not fixed up	Fix up the plug well.
	Switch is out of position	Check the Switch

IX. Jugegment of Fault

Power: 25W
Using time:15m/time

VII. Technological Parameter

There are around 650 points in our body from neck to waist.And every point is connected with some organs or body condition.Massaging can promote blood circulation,promote metabolism and build up the ability agaisit disease.Massager will make you have massaging delight and healthy pleasure by easy ways.



VII. Chart of body point

1. Put the plug into the socket.
2. Press "ON/OFF" button to start or close.
3. Press "VITESSE" button to increase or decrease the massage strength, four speeds to choose.
4. Press "PROGRAMME" button to choose automatic mode,three modes (RELAXANT,ENERGISANT, STIMULANT) to choose.
5. Press "HEAT" button to start or close heating function, the indicator light will be on or off. (Need to fit on heating massage head before open heating function,otherwise no heating effect.)
6. All applications will be automatically closed after 15 minutes.
7. Please pull out the plug after using.

VIIHow to use

Note: Only with heating massage head can be used for heating function



The massager hammer is equiped with three sets massage head accessories.

V. Massage Head Accessories